

# Drain-Guard.


## Percutaneous Catheter Securement Device

Drain-Guard is an anchoring device for percutaneous drainage catheters that are protruding from the patient's skin.



### Product benefits:

- **Reduces risk of accidental displacement of catheter**
- **Safe anchoring and protection of catheter**
- **Prevents contamination of puncture site from environment**
- **Allows inspection of puncture site**
- **Reduces catheter movement**
- **Can prevent kinking of catheter**



### Product features:

Smooth non-woven cover with good adhesion for secure fixation of catheter

Perforated cover for easy replacement of anchoring device

PE-foam support

Thin, comfortable hydrocolloid carrier

Transparent hydrocolloid cover allows monitoring of puncture site

## Safe anchoring of catheters

- Anchoring device for percutaneous drainage catheters that are protruding from the patient's skin in an angle approximately perpendicular to the patient's skin surface
- The device consists of a skin-side surface with a skin friendly adhesive combined with a specially designed foam block that safely anchors catheters from 5 to 16 French.
- Design prevents catheter movement and migration which can enlarge and potentially cause leakage at the exit site
- Soft foam support block can prevent kinking of the catheter.
- Transparent hydrocolloid cover allows continuous monitoring of the puncture site and provides efficient protection against contamination.

## Patient Comfort

- Soft and pliable adhesive carrier offering substantially improved comfort to the patient and reduced need for suturing
- Hydrocolloid has been carefully tested for skin irritation, thus minimising the risk of allergic reaction even when a number of Drain-Guard devices are used consecutively on the same site.
- Device design enables greater freedom of movement for the patient

## Ease of Use

- The Drain-Guard catheter anchoring device is easy to apply to the patient's skin and changing of the device may even be carried out by a family member.

## Ordering Information

Product code	Description	Design type	Product code	Box quantity	PIP Code
AD-8150	Drain-Guard fixation device	Oval - Anti-Kink - Hydrocolloid	AD-8150	25	419-5350

Manufactured by

**Copenhagen  
MedLab**

Distributed by:

*It's Interventional*  
Albreda House  
Lydgate Lane  
Sheffield, S10 5FH

T: +44 (0) 114 268 8880  
F: +44 (0) 114 268 8881  
orders@itsinterventional.com  
www.itsinterventional.com



***It's Interventional.***

*Patient centred people*