

# Epi-Guard.


## Fixation Device for Epidural Catheters

Epi-Guard is a fixation device for epidural catheters. The device consists of a skin-side surface with a skin friendly adhesive combined with a special design.





### Product benefits:

- **Reduces risk of accidental removal of catheter**
- **Safe anchoring and protection of catheter**
- **Allows inspection of puncture site and catheter position**
- **Can prevent kinking of catheter**
- **Can prevent contamination of puncture site from environment**
- **Easy replacement of anchoring device**



### Product features:



## Safe anchoring of catheters

- Anchoring device for epidural catheters
- The device consists of a skin-side surface with a skin friendly adhesive combined with a specially designed foam block that safely anchors epidural catheters
- Hydrocolloid dressing absorbs bodily fluids or anaesthetic given prior to insertion of the epidural catheter
- Design prevents catheter movement and migration which can enlarge and potentially cause leakage at the exit site
- Soft foam support block can prevent kinking of the catheter.

## Patient Comfort

- Epidural catheters are typically secured inadequately by PU-film or other adhesive dressings and accidental displacement of the catheter is frequently seen. Not only will this cause the patient to experience a period without sufficient anaesthesia, but a very uncomfortable insertion of a new catheter is also required. When Epi-Guard is being used, accidental displacement of the epidural catheter will be rare.
- The hydrocolloid carrier has been carefully tested for skin irritation, thus minimising the risk of allergic reaction
- The Epi-Guard gives the patient great freedom to move around without worrying about the catheter.

## Ease of use

- The Epi-Guard is easy and fast to apply to the patient's skin and results in a standardised and safe anchoring of the epidural catheter.

## Ordering Information

Description	Design type	Product code	Box quantity
Epi-Guard fixation device	Oval - Anti-Kink - Hydrocolloid	AD-8170	25

Manufactured by

**Copenhagen  
MedLab**

Distributed by:

*It's Interventional*  
Albreda House  
Lydgate Lane  
Sheffield, S10 5FH

T: +44 (0) 114 268 8880  
F: +44 (0) 114 268 8881  
orders@itsinterventional.com  
www.itsinterventional.com



***It's Interventional.***

*Patient centred people*