

Drainage Bags.



Our range of drainage bags have been developed specifically for patients requiring percutaneous drainage. The luer-lock hub provides a secure connection between the catheter and the drainage bag, maintaining a closed system and reducing chance of fluid leakage or infection. A selection of sizes from 800ml to 2L are available.

800ml Day bag product benefits:

- **Secure connection to the catheter, maintaining a closed system to reduce chance of leakage or infection**
- **Prevents friction & sweating when sat against patients skin**
- **Can be secured in a way that is most comfortable and convenient for the patient, either around the leg or waist**
- **Patient confidence for emptying bag at 1/3 to 1/2 full**
- **Prevents chance of mechanical occlusion and further complications**



Product features:

Hand washing & glove safety signs

Soft fabric backing

Elastic cloth belt

Luer-lock connector

Anti-kink tubing

100ml markings



2L Night bag product benefits:

- **Secure connection to the catheter, maintaining a closed system to reduce chance of leakage or infection**
- **Convenient & easy to use**
- **Easy visualisation for when drainage bag requires emptying**
- **User confidence knowing if tap is open or closed to prevent chance of accidental leakage**



Product features:

Luer-lock connector


Patient information label

Easy to use T-Tap

Hanging threads to attach to night stand

Single use

200ml markings



Ordering Information

Description	Design type	Product code	Box quantity
Drainage bag 800ml	Inc. ML Luer Lock Drainage Line	LU800-DF	25
Drainage bag 2L	Inc. ML Luer Lock Drainage Line	UB-121T-1(ML)-DF	25
Leg bag with fabric backing and Velcro strap 800ml	Inc. ML Luer Lock Drainage Line	DB-800C-MS-V-DF	10

Manufactured by

BIOTEQ

Distributed by:

It's Interventional
Albreda House
Lydgate Lane
Sheffield, S10 5FH

T: +44 (0) 114 268 8880
F: +44 (0) 114 268 8881
orders@itsinterventional.com
www.itsinterventional.com

2022-01-28-Drainage Bags Brochure Sheet



It's Interventional.

Patient centred people