

***It's Interventional.***

*Patient centred people*

# ***MedOne automatic biopsy device.***



**BIOPSY**

*hello@itsinterventional.com*

*itsinterventional.com*



# *It's about accessibility...*

*The MedOne biopsy device features a wingless charging design, paired with an offset needle position to facilitate access to soft tissue lesions. Through optimising spring compression, MedOne delivers good quality, reliable tissue capture even in dense and calcified tissue. Licensed for use in breast, prostate, kidney, lung and liver, the MedOne offers a solution to most soft tissue biopsy needs.*

## **Product benefits:**






- ***"Wingless" handle with an offset needle allows for optimal access.***
- ***Echogenic treatment provides visibility under U/S for all targeted biopsies.***
- ***Side & end fire buttons to cater for clinician choice and procedure requirements.***
- ***Smooth charging mechanism.***
- ***Reliable firing.***
- ***Lightweight compact design.***
- ***Suitable for stereotactic breast biopsies.***
- ***Compatible with prostate specific access systems.***



## **Additional information?**

If you would like further assistance with using MedOne, or need advice on how to order, please call us on: **0114 268 8880**

**Product features:**



## Ordering Information

Product Code	Description	Gauge	Colour	Box Quantity
M0-12-080-00	MedOne 12g x 8cm	12		10
M0-12-100-00	MedOne 12g x 10cm			10
M0-12-120-00	MedOne 12g x 12cm			10
M0-12-160-00	MedOne 12g x 16cm			10
M0-12-180-00	MedOne 12g x 18cm			10
M0-12-200-00	MedOne 12g x 20cm			10
M0-12-250-00	MedOne 12g x 25cm			10
M0-12-300-00	MedOne 12g x 300cm			10
M0-14-080-00	MedOne 14g x 8cm	14		10
M0-14-100-00	MedOne 14g x 10cm			10
M0-14-120-00	MedOne 14g x 12cm			10
M0-14-160-00	MedOne 14g x 16cm			10
M0-14-180-00	MedOne 14g x 18cm			10
M0-14-200-00	MedOne 14g x 20cm			10
M0-14-250-00	MedOne 14g x 25cm			10
M0-14-300-00	MedOne 14g x 30cm			10
M0-16-080-00	MedOne 16g x 8cm	16		10
M0-16-100-00	MedOne 16g x 10cm			10
M0-16-120-00	MedOne 16g x 12cm			10
M0-16-160-00	MedOne 16g x 16cm			10
M0-16-180-00	MedOne 16g x 18cm			10
M0-16-200-00	MedOne 16g x 20cm			10
M0-16-250-00	MedOne 16g x 25cm			10
M0-16-300-00	MedOne 16g x 30cm			10
M0-18-080-00	MedOne 18g x 8cm	18		10
M0-18-100-00	MedOne 18g x 10cm			10
M0-18-120-00	MedOne 18g x 12cm			10
M0-18-160-00	MedOne 18g x 16cm			10
M0-18-180-00	MedOne 18g x 18cm			10
M0-18-200-00	MedOne 18g x 20cm			10
M0-18-250-00	MedOne 18g x 25cm			10
M0-18-300-00	MedOne 18g x 30cm			10
M0-20-080-00	MedOne 20g x 8cm	20		10
M0-20-100-00	MedOne 20g x 10cm			10
M0-20-120-00	MedOne 20g x 12cm			10
M0-20-160-00	MedOne 20g x 16cm			10
M0-20-180-00	MedOne 20g x 18cm			10
M0-20-200-00	MedOne 20g x 20cm			10
M0-20-250-00	MedOne 20g x 25cm			10
M0-20-300-00	MedOne 20g x 30cm			10

Not all codes are stocked items. Non-stock items available on request, lead times may apply. Please contact us for size availability.

Manufactured by



Distributed by:

It's Interventional  
Albreda House  
Lydgate Lane  
Sheffield, S10 5FH

T: +44 (0) 114 268 8880  
F: +44 (0) 114 268 8881  
orders@itsinterventional.com  
www.itsinterventional.com

2023-10-30-MedOne Brochure



**It's Interventional.**

Patient centred people